



**NOBAMED Paul Danz AG**

## **SCHLITZKOMPRESSEN-steril**

**REF 751207**

### **Product Description, Intended use, Application**

Gauze slit swab, **12-ply**, made of 17 thread dressing gauze, **size 7.5 cm x 7.5 cm**, individually **packed per 2 pieces** in sterile packaging, with Y-cut. The gauze slit swab serves as a wound cover of inlets or outlets of probes, drainages, catheters, Kirschner's wire tractions or external fixators. The gauze slit swab is for single use only.

### **Composition**

Cotton

### **Incident reporting**

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

### **Normative and Regulative Requirements, Common Standards**

Medical Device according to *MDD 93/42/EEC*, *MDR (EU) 2017/745*.

The gauze slit swabs comply with the requirements of the DIN EN 14079.

Sterilization of the product complies with DIN EN ISO 11135.

The product does not contain dangerous toxic substances according to REACH.

### **Packaging**

Primary packaging:	paper-paper-packaging
Secondary packaging:	folding box made of cellulose
Tertiary packaging:	carton made of cellulose

### **Symbols used in labelling**

Explanation at [www.nobamed.com](http://www.nobamed.com)



Marking on all packaging levels with CE and according to DIN EN ISO 15223-1- and DIN EN ISO 20417

### **Storage and Transport**

Dry and dustfree, protect from sunlight.

### **Sterile device**

Before using a sterile product, visually inspect the packaging to ensure that it is intact.

### **Single use device**

Reusing a single use medical device can lead to microbiological danger. Reprocessing for reuse can decrease the product's performance significantly.

### **Disposal**

According to locally applicable legal regulations and standards of infection prophylaxis.