



**NOBAMED Paul Danz AG**

## **NOBADRAPE®**

### **Dental drape**

**REF 762750**

#### **Product Description, Intended use, Application**

The 2-ply blue dental drape with diamond-shaped opening (10 x 14 cm) is individually sterile packed. It is used as sterile cover during jaw surgery. The top layer consists of an absorbent polypropylene nonwoven. The bottom layer consists of a moisture-impermeable polyethylene film providing a reliable barrier protection.

Easily removable control stickers on the packaging facilitate surgical documentation.

#### **Composition**

Polypropylene, polyethylene

#### **Contraindications**

The product should not be used in the case of a known allergy against the material.

#### **Normative and Regulative Requirements, Common Standards**

Medical Device according to *MDD 93/42/EEC, MDR (EU) 2017/745*.

The product complies with the requirements of the **DIN EN 13795-1**, "Surgical clothing and drapes – Requirements and test methods – Part 1: Surgical drapes and gowns (requirements of the performance level 'high' and 'standard')".

Sterilization of the product complies with DIN EN 11135.

The product does not contain dangerous toxic substances according to REACH.

#### **Packaging**

Primary:	paper-film packaging
Secondary:	folding box made of cellulose
Tertiary:	carton made of cellulose

#### **Symbols used in labelling**

Explanation at [www.nobamed.com](http://www.nobamed.com)



Marking on all packaging levels with CE and according to DIN EN ISO 15223-1- and DIN EN ISO 20417

#### **Storage and transport**

To be stored in a dry and dust-free environment, protect from sunlight.

#### **Sterile device**

Before using a sterile product, visually inspect the packaging to ensure that it is intact.

#### **Single use device**

Reusing a single use medical device can lead to microbiological danger. Reprocessing for reuse can decrease the product's performance significantly.

#### **Disposal**

According to locally applicable legal regulations and standards of infection prophylaxis.