



NOBAMED Paul Danz AG

VERBANDTUCH-steril

REF 776406

Product Description; Intended use, Application

The wound sheet, **size BR** (small) measuring **40 cm x 60 cm**, is used for aseptic wound coverage and for pressure-free coverage of large-scale injuries. VERBANDTUCH-steril is for single use, especially in first-aid, and is individually sealed in a paper-film packaging.

Composition

Polypropylene

Contra-Indications

The product should not be used in the case of a known allergy against the material.

Notes

DIN 13152-BR 400 x 600 mm

Incident reporting

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Normative and Regulatory Requirements, Common Standards

Medical Device according to *MDD 93/42/EEC, MDR (EU) 2017/745*.

The wound sheet complies with the normative requirements of DIN 13152.

Sterilization in compliance with DIN EN 11135.

The product does not contain dangerous toxic substances according to REACH.

Packaging

Primary packaging: paper-film packaging

Secondary packaging: PE bag

Tertiary packaging: carton made of cellulose

Symbols used in labelling

Explanation at www.nobamed.com



Marking on all packaging levels with CE and according to DIN EN ISO 15223-1- and DIN EN ISO 20417

Storage and Transport

Dry and dustfree, protect from sunlight.

Sterile device

Before using a sterile product, visually inspect the packaging to ensure that it is intact.

Single use device

Reusing a single use medical device can lead to microbiological danger. Reprocessing for reuse can decrease the product's performance significantly.

Disposal

According to locally applicable legal regulations and standards of infection prophylaxis.