



**NOBAMED Paul Danz AG**

## VERBANDTUCH-steril

REF 776608

### Product Description; Intended use, Application

The wound sheet, **size A** (medium) measuring **60 cm x 80 cm**, is used for aseptic wound coverage and for pressure-free coverage of large-scale injuries. VERBANDTUCH-steril is for single use, especially in first-aid, and is individually sealed in a paper-film packaging.

### Composition

Polypropylene

### Contra-Indications

The product should not be used in the case of a known allergy against the material.

### Notes

DIN 13152-A 600 x 800 mm

### Incident reporting

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

### Normative and Regulatory Requirements, Common Standards

Medical Device according to *MDD 93/42/EEC, MDR (EU) 2017/745*.

The wound sheet complies with the normative requirements of DIN 13152.

Sterilization in compliance with DIN EN 11135.

The product does not contain dangerous toxic substances according to REACH.

### Packaging

Primary packaging: paper-film packaging

Secondary packaging: PE bag

Tertiary packaging: carton made of cellulose

### Symbols used in labelling

Explanation at [www.nobamed.com](http://www.nobamed.com)



Marking on all packaging levels with CE and according to DIN EN ISO 15223-1- and DIN EN ISO 20417

### Storage and Transport

Dry and dustfree, protect from sunlight.

### Sterile device

Before using a sterile product, visually inspect the packaging to ensure that it is intact.

### Single use device

Reusing a single use medical device can lead to microbiological danger. Reprocessing for reuse can decrease the product's performance significantly.

### Disposal

According to locally applicable legal regulations and standards of infection prophylaxis.