



NOBAMED Paul Danz AG

NOBALUMENAL®-steril®

REF 873525

Product Description, Intended Use, Application

The sterile eye pad, **size 54 x 70 mm** consists of a soft non-woven and is provided with an absorbent filling. It is closed all around so that no fluffs can emerge. Due to its **oval shape** the pad adapts to the eye socket. The eye pad supports eye dressing as absorbent, padding, protecting or caring part. NOBALUMENAL®-steril is for single use only and individually packed.

Composition

Polyester, viscose, rayon

Contraindications

The product should not be used in the case of a known allergy against the material.

Incident reporting

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Normative and Regulatory Requirements, Common Standards

Medical Device according to *MDD 93/42/EEC*, *MDR (EU) 2017/745*.

Sterilization of the product complies with DIN EN 11135.

The product does not contain dangerous toxic substances according to REACH.

Packaging

Primary packaging:	paper-film packaging
Secondary packaging:	folding box made of cellulose
Tertiary packaging:	carton made of cellulose

Symbols used in labelling

Explanation at www.nobamed.com



Marking on all packaging levels with CE and according to DIN EN ISO 15223-1- and DIN EN ISO 20417

Storage and Transport

Dry and dustfree

Sterile device

Before using a sterile product, visually inspect the packaging to ensure that it is intact.

Single use device

Reusing a single use medical device can lead to microbiological danger. Reprocessing for reuse can decrease the product's performance significantly.

Disposal

According to locally applicable legal regulations and standards of infection prophylaxis.