



NOBAMED Paul Danz AG

NOBADRAINAL[®]-steril

REF 874355

Product Description, Intended Use, Application

NOBADRAINAL[®]-steril is a swab with Y-cut, size **4 cm x 4 cm** and provided with a punched hole. It can therefore be used for hole diameter of 5-10 mm and for probes, catheters, drainages and extensions with an outside diameter of 9–30 Charrière (Ch). The sterile, aluminium-coated, and thereby non-adhesive drainage swab is used for the outlet of probes, catheters, drainages and extensions. NOBADRAINAL[®]-sterile is for single use only and individually packed.

Composition

Polyethylen, aluminium, polyester, rayon

Contraindications

The product should not be used in the case of a known allergy against the material.

Incident reporting

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Normative and Regulative Requirements, Common Standards

Medical Device according to *MDD 93/42/EEC*, *MDR (EU) 2017/745*.

Sterilization of the product complies with DIN EN 11135.

The product does not contain dangerous toxic substances according to REACH.

Packaging

Primary packaging: paper-paper packaging

Secondary packaging: folding box made of cellulose

Tertiary packaging: carton made of cellulose

Symbols used in labelling

Explanation at www.nobamed.com



Marking on all packaging levels with CE and according to DIN EN ISO 15223-1- and DIN EN ISO 20417

Storage and Transport

Dry and dustfree

Sterile device

Before using a sterile product, visually inspect the packaging to ensure that it is intact.

Single use device

Reusing a single use medical device can lead to microbiological danger. Reprocessing for reuse can decrease the product's performance significantly.

Disposal

According to locally applicable legal regulations and standards of infection prophylaxis