



NOBAMED Paul Danz AG

NOBAKARP®

REF 900611

Product Description, Intended use, Application

NOBAKARP® is a disposable patient blanket (size: 110 cm x 190 cm) that is used during patient transport and first aid. It consists of polypropylene non-woven and is filled with 6 layers of cellulose. The blanket is sewn all around and has red fixation corners.

Composition

Polypropylene, cellulose

Contraindications

The product should not be used in the case of a known allergy against the material.

Incident reporting

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Normative and Regulatory Requirements, Common Standards

Medical Device according to MDD 93/42/EEC, MDR (EU) 2017/745.

The product does not contain dangerous toxic substances according to REACH.

Packaging

Primary packaging: PP pouch

Secondary packaging: Carton made of cellulose

Symbols used in labelling

Explanation at www.nobamed.com



Marking on all packaging levels with CE and according to DIN EN ISO 15223-1 and DIN EN ISO 20417

Storage and Transport

Dry and dustfree

Single use device

Reusing a single use medical device can lead to microbiological danger. Reprocessing for reuse can decrease the product's performance significantly.

Disposal

According to locally applicable legal regulations and standards of infection prophylaxis.