



**NOBAMED Paul Danz AG**

## **NOBAKARP®**

**REF 900914**

### **Product Description, Intended use, Application**

NOBAKARP® is a disposable patient winter blanket (size: 110 cm x 190 cm) that is used during patient transport and first aid. It consists of polypropylene non-woven and is filled with layers of cellulose. The blanket is sewn all around and has red fixation corners.

### **Composition**

Polypropylene, cellulose

### **Contraindications**

The product should not be used in the case of a known allergy against the material.

### **Incident reporting**

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

### **Normative and Regulative Requirements, Common Standards**

Medical Device according to MDD 93/42/EEC, MDR (EU) 2017/745.

The product does not contain dangerous toxic substances according to REACH.

### **Packaging**

Primary packaging: PP pouch

Secondary packaging: Carton made of cellulose

### **Symbols used in labelling**

Explanation at [www.nobamed.com](http://www.nobamed.com)



Marking on all packaging levels with CE and according to DIN EN ISO 15223-1 and DIN EN ISO 20417

### **Storage and Transport**

Dry and dustfree

### **Single use device**

Reusing a single use medical device can lead to microbiological danger. Reprocessing for reuse can decrease the product's performance significantly.

### **Disposal**

According to locally applicable legal regulations and standards of infection prophylaxis.