



NOBAMED Paul Danz AG

UNTERSUCHUNGSFINGERLING

REF 903460

Product Description and Purpose

The light-coloured, rolled, powder-free examination fingerstall is made of latex and serves for rectal examination. The thickness is $\geq 0,15$ mm at the sheath and $\geq 0,2$ mm at the fingertip

Composition

Latex

No thiurames, toluenes or ureas are used.

Normative and Legal Requirements

UNTERSUCHUNGSFINGERLING (examination fingerstall) is a medical product according to directive MDD 93/42/EEC and the regulation MDR (EU) 2017/745 and is classified as class I, rule 5.

The virus imperviousness is tested in accordance with **ASTM F1671-07**.

The product does not contain dangerous toxic substances according to REACH. It has CE marking and DIN EN ISO 15223-1 labels on all its packaging.

Packaging

Primary packaging: PP bag

Secondary packaging: carton made of cellulose

Storage

To be stored in a dry and dust-free environment, protected against direct sunlight.

The product bears following labels and symbols:

