



NOBAMED Paul Danz AG

Armtragegurt

REF 910035

Product Description, Intended use, Application

The black arm sling support (width: 35 mm) is provided with a double sling and a clasp which is adjustable. By means of the product the forearm can be immobilized for a short time.

Composition

Cotton, stainless steel

Contraindication

The product should not be used in the case of a known allergy against the material.

Incident reporting

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Normative and Regulative Requirements, Common Standards

Medical Device according to *MDD 93/42/EEC, MDR (EU) 2017/745*.

The product does not contain dangerous toxic substances according to REACH.

Packaging

Primary:	PE bag
Secondary:	carton made of cellulose

Symbols used in labelling

Explanations at www.nobamed.com



Marking on all packaging levels with CE and according to DIN EN ISO 15223-1 and DIN EN ISO 20417

Storage and Transport

Dry and dustfree

Single use device

Intended for use on one patient. The re-processing of single use devices for reuse can decrease the product's performance significantly. Reusing a single use device can lead to microbiological danger.

Disposal

According to locally applicable legal regulations and standards of infection prophylaxis.