



NOBAMED Paul Danz AG

NOBARIGID[®]-closed

REF 943350

Product Description, Intended use, Application

NOBARIGID[®]-closed (Size S) is a ready-made aluminium splint with an inside padding on both sides. It is used for the immobilization of the finger. The splint may be applied without the need of additional adaption in most cases. An adequate bandage is recommended to fix the splint.

Composition

Aluminium, polyester-based polyurethan foam

Contraindications

The product should not be used in the case of a known allergy against the material.

Incident reporting

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Normative and Regulative Requirements, Common Standards

Medical Device according to *MDD 93/42/EEC, MDR (EU) 2017/745*.

The product does not contain dangerous toxic substances according to REACH.

Packaging

Primary packaging: PE bag

Secondary packaging: carton made of cellulose

Symbols used in labelling

Explanation at www.nobamed.com



Storage and transportation

To be stored in a dry and dust-free environment

Single use device

Reusing a single use medical device can lead to microbiological danger. Reprocessing for reuse can decrease the product's performance significantly.

Disposal

According to locally applicable legal regulations and standards of infection prophylaxis.