



NOBAMED Paul Danz AG

SUSPENSORIUM

REF 950004

Product Description, Intended use, Application

The suspensory, size 4= XS, is provided with a detachable bag. It is washable and is used for the fixation of the testicles after operations or in the case of oedema.

Composition

Cotton, clip made of iron, buttons made of polypropylene

Contraindications

The product should not be used in the case of a known allergy against the material.

Incident reporting

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Normative and Legal Requirements

Medical Device according to MDD 93/42/EEC, MDR (EU) 2017/745.

The product does not contain dangerous toxic substances according to REACH.

Sizes

Size 4-6	child
Size 8	adolescent and adult
Size 10-12	seniors

Packaging

Primary: folding box made of cellulose

Secondary: carton made of cellulose

Symbols used in labelling

Explanation at www.nobamed.com



Marking on all packaging levels with CE and according to DIN EN ISO 15223-1 and ISO 20417.

Storage and Transport

Dry and dustfree

Single use device

Intended for use on one patient. Reusing a single use medical device can lead to microbiological danger. Reprocessing for reuse can decrease the product's performance significantly.

Disposal

According to locally applicable legal regulations and standards of infection prophylaxis.