



NOBAMED Paul Danz AG

SICHERHEITSNADELN

REF 962751

Product Description, Intended use, Application

The SAFETY PINS (size 1: 6 mm x 34 mm) are made of stainless steel and are used to fix bandages, like triangular bandages for example. They can also be used to fix Clavicula bandages, Gilchrist bandages and other immobilization bandages. The safety pins are sterilizable by steam (134°C/2 bar).

Composition

Nickel-plated steel

Incident reporting

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Normative and Regulatory Requirements, Common Standards

Accessory for a Medical Device according to MDD 93/42/EEC, MDR (EU) 2017/745.

The product does not contain dangerous toxic substances according to REACH.

Packaging

Primary packaging: folding box made of cellulose

Secondary packaging: film wrapping

Tertiary packaging: carton made of cellulose

Symbols used in labelling

Explanation at www.nobamed.com



Marking on all packaging levels with CE and according to DIN EN ISO 15223-1 and DIN EN ISO 20417

Storage and transport

To be stored in a dry and dust-free environment

Disposal

According to locally applicable legal regulations and standards of infection prophylaxis.