



NOBAMED Paul Danz AG

RUDACLUDE®

REF 050130

Product Description and Purpose

The sterile, oval eye plaster (size L) is provided with a black, opaque layer above which there is a white, soft wound contact layer. RUDACLUDE® has a circumferential adhesive surface so that the eye is optimally protected. The eye plaster can, in sterile condition, be used in the OR. It is especially suited after operations when light must be prevented from shining into the eye for a specific time. Also suited for strabismus (squint) to cover one eye.

Contra-Indications

Known allergy to one of the ingredients listed below.

Composition

Polyester, viscose, polyethylene, polyacrylate adhesive

Normative and Legal Requirements

RUDACLUDE® is a medical product according to the MPG, the directive 93/42/EEC and the regulation MDR (EU 2017/745) and is classified as class Is, according to rule 4.

Sterilization of the product complies with DIN EN 11135.

The product does not contain dangerous toxic substances according to REACH.

It has CE marking and DIN EN ISO 15223-1 labels on all its packaging.

Packaging

Primary packaging:	paper-paper packaging
Secondary packaging:	folding box made of cellulose
Tertiary packaging:	carton made of cellulose

The product bears the following symbols and labels:



Contents sterile unless package is damaged