

seca 201

Ergonomic circumference measuring tape



- Retractable measuring tape for precise body measurements.
- Printed scale on side makes results easy-to-read.
- Convenient one-handed operation - fasten start of tape to retraction mechanism.
- Compact and easy to store.

seca 201:

A high-quality
medical product.

Today's circumference tapes also have to satisfy high medical demands. That is why seca has developed the seca 201 which can measure girth with millimeter precision. The high-quality mechanics guarantee that the 205 cm long tape unwinds easily and locks in place precisely. The housing sits comfortably and securely in the user's hand and is tough enough to withstand accidental dropping.

Retraction mechanism for convenient handling.

The seca 201 is a true medical product. It's very easy to pull out the tape, which just as easily clicks back into place. Handling is extraordinarily convenient. Made of durable materials, the scale features a tape with a top measurement of 205 centimeters and its own low weight of only 50 grams.



Millimeter graduation for
accurate measurements.



High-quality
materials for precise
functionality.

seca 201

Technical Data

- Measurement range: 0–205 cm / 0–80"
- Division: 1 mm / 1/8"
- Diameter: 70 mm / 2,8"
- Height: 22 mm / 0,9"
- Weight of device: 50 g / 0,2 lbs