



**NOBAMED Paul Danz AG**

## NOBAPRON® ROLLE

**REF 902140**

### **Product Description, Intended use, Application**

The white disposable apron (size 140 cm x 80 cm) for use in patient care is made of foil. 50 aprons are on a roll.

The use of the apron reduces the spread of germs and contaminations in patient care and medical activities.

### **Composition**

Polyethylene

### **Incident reporting**

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

### **Normative and Regulative Requirements, Common Standards**

Medical Device according to *MDD 93/42/EEC, MDR (EU) 2017/745*.

The product does not contain dangerous toxic substances according to REACH.

### **Packaging**

Primary: PE bag

Secondary: carton made of cellulose

### **Symbols used in labelling**

Explanation at [www.nobamed.com](http://www.nobamed.com)



Marking on all packaging levels with CE and according to DIN EN ISO 15223-1 and DIN EN ISO 20417

### **Storage and transportation**

dry and dust-free

### **Single use device**

Reusing a single use medical device can lead to microbiological danger. Reprocessing for reuse can decrease the product's performance significantly.

### **Disposal**

According to locally applicable legal regulations and standards of infection prophylaxis.