



NOBAMED Paul Danz AG

NOBAPRON® ROLLE

REF 902160

Product Description, Intended use, Application

The white disposable apron (size 160 cm x 80 cm) for use in patient care is made of foil. 50 aprons are on a roll.

The use of the apron reduces the spread of germs and contaminations in patient care and medical activities.

Composition

Polyethylene

Incident reporting

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Normative and Regulative Requirements, Common Standards

Medical Device according to *MDD 93/42/EEC, MDR (EU) 2017/745*.

The product does not contain dangerous toxic substances according to REACH.

Packaging

Primary: PE bag

Secondary: carton made of cellulose

Symbols used in labelling

Explanation at www.nobamed.com



Marking on all packaging levels with CE and according to DIN EN ISO 15223-1 and DIN EN ISO 20417

Storage and transportation

dry and dust-free

Single use device

Reusing a single use medical device can lead to microbiological danger. Reprocessing for reuse can decrease the product's performance significantly.

Disposal

According to locally applicable legal regulations and standards of infection prophylaxis.