



NOBAMED Paul Danz AG

NOBAHEBAN®

REF 453403

Product Description and Purpose

The brown cohesive compression bandage (size 3.75 cm x 4.5 m) has a stretchability of approx. 100 %. It is used for bandages for short-term compression. It is applied over a gauze ball or swab. It is also used as a supporting and relieving bandage

Note:

During application, cohesive bandages should not be stretched too much, as otherwise the arterial blood supply might be disturbed.

Contraindications

Known allergy against a component listed below.

Composition

Polyamide, Polyester, Latex, Elastomere

Normative and Legal Requirements

NOBAHEBAN® is a medical product according to directive MDD 93/42/EEC and the regulation MDR (EU) 2017/745 and is classified as class I, rule 1.

The product does not contain dangerous toxic substances according to REACH. It has CE marking and DIN EN ISO 15223-1 labels on all its packaging.

Packaging

Primary packaging:	foil
Secondary packaging:	folding box made of cellulose
Tertiary packaging:	carton made of cellulose

Storage

To be stored in a dry and dust-free environment.

The product bears the following marking and symbols.

