



NOBAMED Paul Danz AG

APPLIKATOR

REF 320005

Product Description, Intended use, Application

The APPLICATOR is a device made of stainless steel (5 cm, for big finger/ toe), consisting of two rings connected with longitudinal rods. The APPLICATOR is used for easy and quick application of tubular bandages (such as NOBATRICOT®) during wound care.

Composition

Stainless steel

Incident reporting

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Normative and Regulative Requirements, Common Standards

Accessory for medical device tubular bandage according to the MDD 93/42/EEC and MDR (EU) 2017/745.

The product does not contain dangerous toxic substances according to REACH.

Packaging

Primary packaging:	PE-bag
Secondary packaging:	carton made of cellulose

Symbols used in labelling

Explanations at www.nobamed.com



Marking on all packaging levels with CE and according to DIN EN ISO 15223-1 and ISO 20417.

Storage and Transport

To be stored in a dry and dust-free environment

Disposal

According to locally applicable legal regulations and standards of infection prophylaxis.