



NOBAMED Paul Danz AG

PFLEGESCHÜRZE

REF 664500

Product Description and Purpose

The yellow apron for use in patient care is made of nonwoven. It is available in one universal size. The apron is pulled over the head and the waist ties are fixed. Thumb-loops and elastic cuffs facilitate the easy and tight donning of disposable gloves. The use of the apron reduces the spread of germs and contaminations in patient care and medical activities (minimum exposure to liquids). The disposable apron is for single use.

Contra-Indications

No known contra-indications.

The product should not be used if there is a known hypersensitivity to one of the below-mentioned components.

Composition

Polypropylene

Normative and Legal Requirements

PFLEGESCHÜRZE (protective apron for use in patient care) is a class I, rule 1 medical device according to the Directive 93/42/EEC and the Regulation MDR (EU) 2017/745.

The product does not contain dangerous toxic substances according to REACH. It has CE marking and DIN EN ISO 15223-1 labels on all its packaging.

Packaging

Primary packaging: PE bag

Secondary packaging: carton made of cellulose

Storage

To be stored in a dry and dust-free environment

The product bears following labels and symbols:

