



NOBAMED Paul Danz AG

RUDAVAL[®]

REF 052120

Product Description and Purpose

Sterile, oval, individually packed eye plaster made of a soft nonwoven carrier material and a central wound pad. The circumferential adhesive surface of the plaster provides the necessary protection to the eye.

Contraindications

Known allergy to one of the ingredients listed below.

Composition

Polyester, viscose, polyethylene, polyacrylate adhesive

Normative and Legal Requirements

RUDAVAL[®] is a medical product according to the MPG, the directive 93/42/EEC and the regulation MDR (EU 2017/745) and is classified as class Is, according to rule 4.

Sterilization of the product complies with DIN EN 11135.

The product does not contain dangerous toxic substances according to REACH.

It has CE marking and DIN EN ISO 15223-1 labels on all its packaging.

Packaging

Primary packaging:	paper-paper packaging
Secondary packaging:	folding box made of cellulose
Tertiary packaging:	carton made of cellulose

The product bears the following symbols and labels:

