

metpak

pressure infusion device



Pressure infusion devices are an essential item in any hospital! All the more so when it comes to the high-quality products from **Riester**, which excel by their particularly high degree of **reliability and durability**. Pressure cuffs in combination with a manometer for producing pressure provide **for optimal pressure infusions**. For solutions or blood bags from 500 ml, 1000 ml, 3000 ml or 5000 ml.

- › Sturdy cotton cuff in blue, with suspension loop for infusion stands, washable up to 60°C
- › Precise reading of the level in the bag thanks to a transparent net on the cuff
- › Chrome-plated manometer with easy-to-read 49 mm diameter aluminium scale up to 300 mm Hg
- › Specially hardened copper-beryllium membrane
- › No zero point setting required
- › Latex-free pad, resistant to pressure
- › Latex-free inflation bulb

| | |
|---------------------|---------------------------|
| No. 5270 | metpak 500ml, latex-free |
| No. 5275 | metpak 1000ml, latex-free |
| No. 5270-536 | metpak 3000ml, latex-free |
| No. 5270-537 | metpak 5000ml, latex-free |

Please visit [riester.de](https://www.riester.de) for more information