



NOBAMED Paul Danz AG

HAUBE CLIP CAP

REF 674901

Product Description and Purpose

The caps are made of a light, white spun-bonded nonwoven. The material is permeable to air and therefore breathable. This prevents heat accumulation, which increases the wearing comfort.

The clip cap is used to cover the hair of staff, patients and visitors. A contamination by hair is avoided.

Composition

Polypropylene, Elastomer thread (latex-free)

Normative and Legal Requirements

The Clip Cap is a medical product according to the MPG, the directive 93/42/EEC and the regulation MDR (EU 2017/745) and is classified as class I, according to rule 1.

The product does not contain dangerous toxic substances according to REACH.

It has CE marking and DIN EN ISO 15223-1 labels on all its packaging.

Packaging

| | |
|----------------------|-------------------------------|
| Primary packaging: | folding box made of Cellulose |
| Secondary packaging: | carton made of cellulose |

Storage

To be stored in a dry and dust-free environment

The product bears following labels and symbols:

