



Product Specification – Focus - Saline Eye Wash - 20ml

Product Information

REF:	6058	
Description:	Focus - Saline Eye Wash - 20ml	
Manufacturer:	Longbow First Aid Products Manufactory	
Address:	2/F, Area C, Hantian Industrial Park, Guiping Road Guicheng Subdistrict Nanhai District, 528200 Foshan City, Guangdong Province, China	
EC Rep:	Shanghai International Holding Corp Eiffestrasse 80, 20537, Hamburg, Germany.	
EU Importer:	Reliance Medical Ireland Unit 17, Westlink Industrial Estate, Kylemore Road, Dublin 10, Ireland	
Country of Origin:	China	
Commodity Code:	3004900000	
Shelf Life:	5 years	
Storage Conditions:	Store in a cool, dry place. Store between 5-35°C. Keep out of direct sunlight. Keep Dry.	

Intended Use:	For minor eye irritations.
Medical Device Classification:	Is
Compliance Markings:	CE, Medical Device
Main Materials:	Sodium Chloride 0.9%, Aqua
Product Dimensions:	11.5cmH x 1.5cmW x 2cmL
Product Net Weight:	28.2g
Sterilisation method:	Sterile R
Latex Information:	Latex-free
Number of Uses:	Single application
Warnings or Contra-Indications:	Do not use if package is open or damaged
Disposal Information:	No special disposal method required

The information contained in this product specification is accurate for the above revision date.

However, components are subject to change without notice, when another specification will be issued.



Packaging Information

Packaging Style	SINGLE	PACK/BOX	CASE
Packaging Material(s)	PE 5.5g	Card 28.5g	Card 400g
Recyclable Packaging	Yes	Yes	yes
UN Code(s)	N/A		
Hazard Label(s)	N/A		

Shipping Information

UNIT	QUANTITY	BARCODE	PACKAGING DIMENSIONS			GROSS WEIGHT
			mmL	mmW	mmH	
Single Unit:	Single	N/A	115	20	15	28.5g
Pack/Box:	25	5032483001569	119	90	11	700g
Case:	10 boxes of 25	15032483001566	370	245	140	6.4kg

*The information contained in this product specification is accurate for the above revision date.
However, components are subject to change without notice, when another specification will be issued.*